

# Omne Trinum Perfectum 1.0

---

- Gravity **12.4 BLG**
- ABV **5 %**
- IBU **30**
- SRM **7**
- Style **American Pale Ale**

## Batch size

- Expected quantity of finished beer **20 liter(s)**
- Trub loss **5 %**
- Size with trub loss **23 liter(s)**
- Boil time **60 min**
- Evaporation rate **12 %/h**
- Boil size **28 liter(s)**

## Mash information

- Mash efficiency **80 %**
- Liquor-to-grist ratio **3.5 liter(s) / kg**
- Mash size **17.5 liter(s)**
- Total mash volume **22.5 liter(s)**

## Steps

- Temp **68 C**, Time **60 min**
- Temp **72 C**, Time **15 min**
- Temp **80 C**, Time **1 min**

## Mash step by step

- Heat up **17.5 liter(s)** of strike water to **74.9C**
- Add grains
- Keep mash **60 min** at **68C**
- Keep mash **15 min** at **72C**
- Keep mash **1 min** at **80C**
- Sparge using **15.5 liter(s)** of **76C** water or to achieve **28 liter(s)** of wort

## Fermentables

Type	Name	Amount	Yield	EBC
Grain	Weyermann - Pale Ale	3 kg (60%)	80 %	7
Grain	Weyermann - Monachijski I	1 kg (20%)	80 %	14
Grain	Weyermann - Caraamber	0.5 kg (10%)	80 %	65
Grain	Weyermann - Zakwaszający	0.25 kg (5%)	80 %	5
Grain	Płatki jęczmienne	0.25 kg (5%)	80 %	1

## Hops

Use for	Name	Amount	Time	Alpha acid
Dry Hop	Mosaic (USA - 2018) - burzliwa w hopsiderze	25 g	3 day(s)	13.2 %
Dry Hop	Citra (USA - 2018) - burzliwa w hopsiderze	25 g	3 day(s)	12.5 %
Dry Hop	Cascade (USA - 2018) - burzliwa w hopsiderze	25 g	3 day(s)	5.8 %

Boil	Mosaic (USA - 2018)	25 g	5 min	13.2 %
Boil	Citra (USA - 2018)	25 g	10 min	12.5 %
Boil	Cascade (USA - 2018)	25 g	30 min	5.8 %

## Yeasts

Name	Type	Form	Amount	Laboratory
Safale US-05	Ale	Dry	11.5 g	Fermentis

Starter 1L wody + 80g jasnego ekstraktu słodowego suchego Bruntal na 24h. Fermentacja w 18C.

## Extras

Type	Name	Amount	Use for	Time
Water Agent	Gips piwowarski	4 g	Mash	76 min
Dodane razem z ziarnem.				
Water Agent	Kwas fosforowy (woda do wysładzania)	5 g	Mash	76 min
ml nie gramy.				
Fining	Whirlfolc T	2.5 g	Boil	5 min
Jedna tabletką.				
Other	Askorbinian sodu	5 g	Bottling	---
Rozpuszczone w 100ml zimnej wcześniej przegotowanej wody.				